



INDIRA GANDHI UNIVERSITY, MEERPUR, REWARI

(A State University Established under Haryana Act No. 29 of 2013)

No. IGU/Zoology/2020/138-47

Dated: 27.07.2020

To

Subject: - Sealed quotation for PCR with 1 KVA Online UPS Due on 12.08.2020

Dear Sir/Madam,

Please quote for the article mentioned below so as to reach the undersigned by 12.08.2020 (04:00 PM). The quality, make, packaging and other particulars of each item should be stated.

WHILE SUBMITTING THE QUOTATIONS PLEASE NOTE THE POINTS MENTIONED BELOW:-

1. The quotations are to be sent in a sealed envelope marked QUOTATION DUE ON 12-08-2020 and this envelope may be put in another envelope addressed to the undersigned.
2. The rates of S.T/C.S.T./GST, EXCISE, CUSTOM and another taxes chargeable must be specified.
3. The rates of Insurance, if any, should be specified. The firm will be required to submit original Payee Receipt along with the bill.
4. Please state the time period within which the items will be supplied.
5. GUARANTEE PERIOD OF THE INSTRUMENT MUST BE MENTIONED.
6. Discount if any, is to be given with the quotation
7. The rates quoted should be FOR destination IGU, Meerpur, Rewari.
8. Please state, if you are the original manufacturer/sole distributor/ Authorized dealer of the items quoted by you.

Sr. No.	Name of item	Quantity	Make	Model Cat. No.	Other Specification
1.	PCR with 1 KVA UPS	One			<ol style="list-style-type: none">1. Memory-500 typical programs on board, unlimited with USB flash drive expansion2. Temperature control modes- Calculated and block3. PCR license- Yes4. Programming options- Step based graphical5. Reporting- Exportable run logs, system logs6. Instant incubation- Yes7. Capacity-96 x 0.2 ml tubes, 0.2 ml tube strips, or 1 x 96- well plate8. Maximum ramp temperature-4 °C/sec9. Average ramp temperature-2.5 °C/sec10. Temperature range-2.5 °C/sec11. Temperature accuracy-0.5 °C/sec12. Temperature uniformity-0.5 °C well to well (within row) within 30 sec13. Thermal gradient14. Gradient capacity- Yes, All blocks gradient15. Gradient accuracy-0.5 °C/sec of programmed temperature16. Row uniformity-0.5 °C well to well (within row) within 30 sec17. Gradient range-30-100 °C18. Temperature differential range-1-25 °C19. Dynamic Ramping technology- Adjust the cyler ramp rate, ensuring that the incubation times across the eight rows during a gradient step are identical20. Consumables- Compatible with standard full height tubes, tube strips and 96 well plates

Thanking you,

Yours sincerely,

Chairperson

CC to:-

A copy of the above is forwarded to incharge University website for uploading the same on University website.

1. **M/S Gagan Scientific Centre**
Near shivaji Colony Chowk, I
Inside flyover, Rohtak 124001

2. **M/s Unique Scientific Agency**
Shop No-7, Munjal Complex
Near Ashoka Choke, Rohtak-124001

3. **Biocraft Scientific System Ltd**
Unit 1:37/277 Naglapadi
New Agra -282005

4. **Haryana Scientific & Engineering Corporation**
Rohtak, Sonipat road,
Vikas nagar, Rohtak 124001

5. **M/S Laboratory Instruments and Chemicals**
Cross Road No. 5, Near Classic Hotel
Grain Market, Ambala Cantt-133001

6. **M/S Satguru Enterprise**
3, Vishwakarma Market
Behind Elite Cinema
Hisar-125004

7. **M/S Galaxy Lab. Instruments**
2157 T-1 Guru Arjun Nagar first floor opposite Sahani Hotel
West patel Nagar, New Delhi 110008

8. **M/S The Bharat Instrument and Chemicals**
Shop No-3, Hnd floor, Top Edge Plaza City
Thana Road, Hisar 125001

9. **M/S A.S. Scientific Instrument and Chemical**
Shop No-1496/4 Ward -31, Kamla Nagar
Rohtak-124001

10. **Proactive Healthcare Pvt. Ltd,**
16/6 Lower Ground Floor, Kalka Ji,
New Delhi, Phone-9811544334
Pin Code-110019

g. b. m.
27.7.2020